

Functional Area of Developmental Molecular Biology

Centrifugation and Ultracentrifugation



The **Centrifugation and Ultracentrifugation** platform is equipped with various instruments which allow centrifugation at different speeds (up to 90000 rpm in the case of Ultracentrifugation), with fixed angle or swinging bucket rotors, at controlled temperature.

Manager: *Alessandro Amoroso*
alessandro.amoroso@szn.it

Second Floor - EAST Wing, Landing

Bookable instruments	Description	
Ultracentrifuge BECKMAN Optima XL 100K	Refrigerated centrifuge, speed up to 90.000 rpm. Rotors available: Ti 90- Ti75- Ti70- Ti65.5- Ti60- SW60Ti- SW41Ti- SW25.1- SW28	
Centrifuge BECKMAN Avanti J25	Refrigerated centrifuge for speed up to 25.000 rpm. Rotors available: JA25.50- JS7.5 swing-out	
Centrifuge BECKMAN J2 MC	Refrigerated Centrifuge for speed up to 25000 rpm. Rotors available: JA10- JA14- JA20.	